
amtech®

cart2count

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Battery-powered LAF trolley
with continuous monitoring

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AM Instruments offers systems, facilities and equipment for cleanrooms. AMTech manufactures standard and custom-made products through a careful analysis of the customer's requests and special needs. Our team is able to offer solutions even for unusual and complex dimensions, functions and layouts, thanks to the support of the highly modern design and production processes.

FROM THE DESIGN STAGE TO THE PROJECT STAGE AND THROUGH TO COMPLETION, AM INSTRUMENTS IS PROFESSIONAL, COMPETENT AND PROACTIVE.

With its own multifunctional laboratory featuring dedicated working areas for the maintenance, repair and calibration of instruments, AM Instruments offers assistance and services which are among the most professional on the market.

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cart2count

THE IDEAL SOLUTION FOR SAFE TRANSFER

Cart2Count has been designed to allow the transfer of materials between two Grade A areas, passing through a Grade B area, without the need for connection to the mains electrical supply and with continuous, on-board particle monitoring.

With **Cart2Count**, AM Tech by AM Instruments has made full use of its experience in building laminar flow systems as well as monitoring systems

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New design with continuous particle counter:

- sealed doors and recirculating airflow
- sampling probe in the "worst case" position
- compliance with data integrity in the management of monitoring data
- complete control of process parameters through the integrated comfort panel
- transmission of Wi-Fi data with Modbus TCP for integration with SCADA supervision systems

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The **LAF trolley** has been developed to allow the completely safe handling and transfer of materials between Grade A areas, passing through a Grade B area, without the need for an electrical supply (battery life of up to one hour).

The LAF trolley is commonly used to transfer partially stoppered vials from the filling line to the freeze-dryers, and then from the freeze-dryers to the vial-capping machines.

Air is drawn from the environment through pre-filters located at the sides, filtered by the HEPA filter, forced horizontally into the working area, then expelled through openings located on the doors.

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Depending on the needs of the client, **Cart2Count** can be fitted with a continuous particle monitoring system.

Lighthouse APEX R5P - Ethernet with built-in pump. The particles are counted through a classic Light Scattering system, and the light source for the counting is a laser diode with a threshold of 0.5 microns for the first channel and 5.0 microns for the second, in the standard model.

Compatible with Hydrogen Peroxide (H₂O₂) vapour in concentrations of up to 1000 ppM. The appliance features buffer storage capable of holding up to 3,000 samples in the event that communication with the data acquisition system is interrupted.

- 0.5 - 5.0 µm standard
- flow rate: 1.0 CFM (28.3 LPM)
- data storage: 3,000 records
- H₂O₂ vapour decontamination compatible
- self-diagnostics for laser, flow and detector
- NIST traceable meets ISO 21501 - 4 and JJF 1190
- view real-time data from your Mac, PC, tablet or smartphone
- communication modes: RS-485, Ethernet, POE, Modbus TCP, ASCII, RTU
- smart-port for configuration, remote display and environmental sensors

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Preparation and implementation of protocols and Risk Analysis reports, FAT, IQ/OQ of the system consisting of LAF Trolley + particle counter + management SCADA.

RISK ANALYSIS

Theoretical analysis with the preparation of a Risk Assessment for the positioning of the isokinetic particle count probe.


FAT

Factory Acceptance Testing at the AM instruments premises for the dimensional and functional verification of the system. Installation Qualification and Operational Qualification.

DOCUMENTATION

Cart2Count is supplied with the following documentation:

- operating manual
- CE certificate of conformity
- original certificates for the HEPA fan and filters
- equipment and instrument list with parts descriptions
- electrical component list with descriptions
- spare parts list with relevant ID numbers
- technical drawings
- electrical circuit diagrams
- drawings of instruments and piping
- schematic diagrams

Case study

DEVIAZIONE DEVIATION

“No continuous monitoring of particles is carried out throughout the critical process for A class areas. Transportation of not finally stoppered vials with products through B class area is carried out in portable untight carts with horizontal class A LAF, no continuous particle monitoring is carried out in the carts. The particles are checked only before transportation beginning and after product unloading.”

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SOLUTION

During the design phase the most critical aspects were examined, for which highly satisfactory solutions were reached: minimising stagnation and turbulence zones in the area where the trays are loaded by means of sealed doors and a recirculating airflow system; carrying out a risk assessment to install the isokinetic probe at the most relevant point; making data transmission “secure” and “compliant” by installing dedicated access points; making changes to the Wi-Fi system to ensure a constant connection that allows for continuity of data monitoring and integrity. In the event of a server connection error, the PLC on the trolley guarantees data storage. Once the connection is re-established, the system automatically downloads the data to the server. In the event of a failure to communicate with SCADA or full memory (during local storage), a special alarm is triggered. The trolley features a control panel used for self-diagnostics, displaying the battery status, and checking the status of sampling and of the particle counter. The trolley was designed and built and all the verification tests necessary for its release were carried out within the short time frame imposed by the inspection request.

aminstruments

AM Instruments srl
Via Isonzo, 1/C
20812 Limbiate (MB)
Tel. +39 02 872892.1

aminstruments.com
info@aminstruments.com

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001:2015 =